

Product datasheet for **GA100116**

Human AP2B1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AP2B1 CRISPRa kit - CRISPR gene activation of human adaptor related protein complex 2 subunit beta 1
Symbol:	AP2B1
Components:	GA100116G1, AP2B1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100116G2, AP2B1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100116G3, AP2B1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P63010</u>
Synonyms:	ADTB2; AP2-BETA; AP105B; CLAPB1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)