

Product datasheet for **GA100130**

Human Angiotensin II Type 1 Receptor (AGTR1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AGTR1 CRISPRa kit - CRISPR gene activation of human angiotensin II receptor type 1
Symbol:	AGTR1
UniProt ID:	P30556
Synonyms:	AG2S; AGTR1B; AT1; AT1AR; AT1B; AT1BR; AT1R; AT2R1; HAT1R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)