

Product datasheet for **GA100187**

Human AMPD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AMPD3 CRISPRa kit - CRISPR gene activation of human adenosine monophosphate deaminase 3
Symbol:	AMPD3
UniProt ID:	Q01432
Synonyms:	adenosine monophosphate deaminase (isoform E); adenosine monophosphate deaminase 3; AMP aminohydrolase; AMP deaminase 3; erythrocyte-specific AMP deaminase; erythrocyte type AMP deaminase; myoadenylate deaminase
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)