

## Product datasheet for **GA100251**

### Human Prostate Specific Antigen (KLK3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	KLK3 CRISPRa kit - CRISPR gene activation of human kallikrein related peptidase 3
Symbol:	KLK3
UniProt ID:	<a href="#">P07288</a>
Synonyms:	APS; hK3; KLK2A1; PSA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)