

Product datasheet for **GA100257**

Human Aquaporin 3 (AQP3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AQP3 CRISPRa kit - CRISPR gene activation of human aquaporin 3 (Gill blood group)
Symbol:	AQP3
UniProt ID:	Q92482
Synonyms:	AQP-3; GIL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)