

Product datasheet for **GA100263**

Human Androgen Receptor (AR) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AR CRISPRa kit - CRISPR gene activation of human androgen receptor
Symbol:	AR
UniProt ID:	P10275
Synonyms:	AIS; AR8; DHTR; HUMARA; HYSPI; KD; NR3C4; SBMA; SMAX1; TFM
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)