

Product datasheet for **GA100265**

Human A RAF (ARAF) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARAF CRISPRa kit - CRISPR gene activation of human A-Raf proto-oncogene, serine/threonine kinase
Symbol:	ARAF
UniProt ID:	P10398
Synonyms:	A-RAF; ARAF1; PKS2; RAFA1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)