

Product datasheet for **GA100268**

Human Amphiregulin (AREG) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AREG CRISPRa kit - CRISPR gene activation of human amphiregulin
Symbol:	AREG
UniProt ID:	P15514
Synonyms:	AR; AREGB; CRDGF; SDGF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)