

Product datasheet for **GA100276**

Human ARG2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARG2 CRISPRa kit - CRISPR gene activation of human arginase 2
Symbol:	ARG2
UniProt ID:	P78540
Synonyms:	arginase, type II; arginase 2; kidney arginase; L-arginine amidinohydrolase; L-arginine ureahydrolase; nonhepatic arginase
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)