

Product datasheet for **GA100283**

Human ARHGAP4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARHGAP4 CRISPRa kit - CRISPR gene activation of human Rho GTPase activating protein 4
Symbol:	ARHGAP4
Components:	GA100283G1, ARHGAP4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100283G2, ARHGAP4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100283G3, ARHGAP4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P98171
Synonyms:	C1; p115; RGC1; RhoGAP4; SrGAP4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

