

Product datasheet for **GA100296**

Human Arrestin C (ARR3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARR3 CRISPRa kit - CRISPR gene activation of human arrestin 3
Symbol:	ARR3
Components:	GA100296G1, Arrestin C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100296G2, Arrestin C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100296G3, Arrestin C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P36575</u>
Synonyms:	ARRX
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

