

Product datasheet for **GA100297**

Human beta Arrestin 1 (ARRB1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARRB1 CRISPRa kit - CRISPR gene activation of human arrestin beta 1
Symbol:	ARRB1
UniProt ID:	P49407
Synonyms:	ARB1; ARR1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

