


Product datasheet for **GA100306**

Human ART3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ART3 CRISPRa kit - CRISPR gene activation of human ADP-ribosyltransferase 3
Symbol:	ART3
Components:	GA100306G1, ART3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100306G2, ART3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100306G3, ART3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q13508
Synonyms:	ARTC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)