

Product datasheet for **GA100307**

Human ART4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ART4 CRISPRa kit - CRISPR gene activation of human ADP-ribosyltransferase 4 (Dombrock blood group)
Symbol:	ART4
Components:	GA100307G1, ART4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100307G2, ART4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100307G3, ART4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q93070</u>
Synonyms:	ARTC4; CD297; DO; DOK1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)