

Product datasheet for **GA100317**

Human GET3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GET3 CRISPRa kit - CRISPR gene activation of human guided entry of tail-anchored proteins factor 3, ATPase
Symbol:	GET3
UniProt ID:	O43681
Synonyms:	ARSA-I; ARSA1; ASNA-I; GET3; hASNA-I; TRC40
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)