

## Product datasheet for **GA100322**

### Human **ASTN1** activation kit by **CRISPRa**

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                            |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| <b>Description:</b>  | ASTN1 CRISPRa kit – CRISPR gene activation of human astrotactin 1                                                                                                                                                                                                                                          |
| <b>Symbol:</b>       | ASTN1                                                                                                                                                                                                                                                                                                      |
| <b>Components:</b>   | GA100322G1, ASTN1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA100322G2, ASTN1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA100322G3, ASTN1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">O14525</a></u>                                                                                                                                                                                                                                                                              |
| <b>Synonyms:</b>     | ASTN                                                                                                                                                                                                                                                                                                       |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:

