

Product datasheet for **GA100344**

Human SERCA3 (ATP2A3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ATP2A3 CRISPRa kit - CRISPR gene activation of human ATPase sarcoplasmic/endoplasmic reticulum Ca ²⁺ transporting 3
Symbol:	ATP2A3
UniProt ID:	Q93084
Synonyms:	SERCA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)