

Product datasheet for **GA100347**

Human ATP2B3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ATP2B3 CRISPRa kit - CRISPR gene activation of human ATPase plasma membrane Ca ²⁺ transporting 3
Symbol:	ATP2B3
UniProt ID:	Q16720
Synonyms:	CFAP39; CLA2; OPCA; PMCA3; PMCA3a; SCAX1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)