

Product datasheet for GA100360

Human ATP5MC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ATP5MC3 CRISPRa kit - CRISPR gene activation of human ATP synthase membrane subunit c locus 3
Symbol:	ATP5MC3
UniProt ID:	P48201
Synonyms:	MGC125738
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)