

Product datasheet for **GA100373**

Human ATP7A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ATP7A CRISPRa kit - CRISPR gene activation of human ATPase copper transporting alpha
Symbol:	ATP7A
UniProt ID:	Q04656
Synonyms:	DSMAX; MK; MNK; SMAX3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)