


Product datasheet for **GA100378**

Human KIF1A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIF1A CRISPRa kit - CRISPR gene activation of human kinesin family member 1A
Symbol:	KIF1A
Components:	GA100378G1, KIF1A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100378G2, KIF1A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100378G3, KIF1A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q12756</u>
Synonyms:	ATSV; C2orf20; HSN2C; MRD9; SPG30; UNC104
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)