

Product datasheet for **GA100384**

Human AVPR V2 (AVPR2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AVPR2 CRISPRa kit - CRISPR gene activation of human arginine vasopressin receptor 2
Symbol:	AVPR2
UniProt ID:	P30518
Synonyms:	ADHR; DI1; DIR; DIR3; NDI; V2R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)