

Product datasheet for **GA100396**

Human **BAI3 (ADGRB3)** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	ADGRB3 CRISPRa kit - CRISPR gene activation of human adhesion G protein-coupled receptor B3
Symbol:	ADGRB3
UniProt ID:	<u>O60242</u>
Synonyms:	BAI3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)