

Product datasheet for **GA100418**

Human BCMA (TNFRSF17) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF17 CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 17
Symbol:	TNFRSF17
UniProt ID:	Q02223
Synonyms:	BCM; BCMA; CD269; TNFRSF13A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)