

## Product datasheet for **GA100419**

### Human HCN2 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	HCN2 CRISPRa kit - CRISPR gene activation of human hyperpolarization activated cyclic nucleotide gated potassium and sodium channel 2
Symbol:	HCN2
UniProt ID:	<a href="#">Q9UL51</a>
Synonyms:	BCNG-2; BCNG2; HAC-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)