

Product datasheet for **GA100428**

Human Filensin (BFSP1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	BFSP1 CRISPRa kit - CRISPR gene activation of human beaded filament structural protein 1
Symbol:	BFSP1
UniProt ID:	<u>Q12934</u>
Synonyms:	CP94; CP115; CTRCT33; LIFL-H
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)