

## Product datasheet for **GA100457**

### Human POLR3D activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	POLR3D CRISPRa kit - CRISPR gene activation of human RNA polymerase III subunit D
Symbol:	POLR3D
UniProt ID:	<u><a href="#">P05423</a></u>
Synonyms:	BN51T; RPC4; RPC53; TSBN51
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)