

## Product datasheet for **GA100467**

### Human BPI activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | BPI CRISPRa kit - CRISPR gene activation of human bactericidal permeability increasing protein |
| Symbol:       | BPI  |
| UniProt ID:   | <a href="#">P17213</a>   |
| Synonyms:     | BPIFD1; rBPI   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                        |

#### Product images:



[View online >](#)