

Product datasheet for **GA100473**

Human TIS11D (ZFP36L2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZFP36L2 CRISPRa kit - CRISPR gene activation of human ZFP36 ring finger protein like 2
Symbol:	ZFP36L2
Components:	GA100473G1, TIS11D gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100473G2, TIS11D gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100473G3, TIS11D gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P47974</u>
Synonyms:	BRF2; ERF-2; ERF2; RNF162C; TIS11D
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

