

Product datasheet for **GA100494**

Human C1QC activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | C1QC CRISPRa kit - CRISPR gene activation of human complement C1q C chain |
| Symbol: | C1QC |
| UniProt ID: | <u>P02747</u> |
| Synonyms: | C1Q-C; C1QG |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)