

Product datasheet for **GA100499**

Human C3AR1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C3AR1 CRISPRa kit - CRISPR gene activation of human complement C3a receptor 1
Symbol:	C3AR1
UniProt ID:	Q16581
Synonyms:	AZ3B; C3AR; HNFAG09
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)