

Product datasheet for **GA100508**

Human Complement C7 (C7) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | C7 CRISPRa kit - CRISPR gene activation of human complement C7 |
| Symbol: | C7 |
| UniProt ID: | P10643 |
| Synonyms: | complement component 7; complement component C7 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)