


## Product datasheet for **GA100518**

### Human MRF (MYRF) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | MYRF CRISPRa kit - CRISPR gene activation of human myelin regulatory factor  |
| Symbol:       | MYRF   |
| Components:   | GA100518G1, MRF gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA100518G2, MRF gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA100518G3, MRF gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9Y2G1</a></u>  |
| Synonyms:     | C11orf9; MRF; Ndt80; pqn-47  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)  |

#### Product images:



[View online »](#)