

Product datasheet for **GA100525**

Human TMEM50B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM50B CRISPRa kit - CRISPR gene activation of human transmembrane protein 50B
Symbol:	TMEM50B
Components:	GA100525G1, TMEM50B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100525G2, TMEM50B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100525G3, TMEM50B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P56557</u>
Synonyms:	C21orf4; HCVP7TP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)