

Product datasheet for **GA100526**

Human MPPED1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MPPED1 CRISPRa kit – CRISPR gene activation of human metallophosphoesterase domain containing 1
Symbol:	MPPED1
Components:	GA100526G1, MPPED1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100526G2, MPPED1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100526G3, MPPED1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q15442</u>
Synonyms:	239AB; C22orf1; FAM1A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

