

## Product datasheet for **GA100538**

### Human CACNA1A activation kit by CRISPRa

#### Product data:

|               |                                                                                                      |
|---------------|------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                         |
| Description:  | CACNA1A CRISPRa kit - CRISPR gene activation of human calcium voltage-gated channel subunit alpha1 A |
| Symbol:       | CACNA1A                                                                                              |
| UniProt ID:   | <u><a href="#">Q00555</a></u>                                                                        |
| Synonyms:     | APCA; BI; CACNL1A4; CAV2.1; EA2; FHM; HPCA; MHP; MHP1; SCA6                                          |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                              |

#### Product images:



[View online >](#)