

Product datasheet for **GA100541**

Human CaV1.3 (CACNA1D) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CACNA1D CRISPRa kit - CRISPR gene activation of human calcium voltage-gated channel subunit alpha1 D
Symbol:	CACNA1D
UniProt ID:	Q01668
Synonyms:	CACH3; CACN4; CACNL1A2; Cav1.3; CCHL1A2; PASNA; SANDD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)