


Product datasheet for **GA100543**

Human CACNA1F activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CACNA1F CRISPRa kit - CRISPR gene activation of human calcium voltage-gated channel subunit alpha1 F
Symbol:	CACNA1F
UniProt ID:	O60840
Synonyms:	AIED; Cav1.4; Cav1.4alpha1; COD3; COD4; CORDX; CORDX3; CSNB2; CSNB2A; CSNBX2; JM8; JMC8; OA2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)