

Product datasheet for **GA100546**

Human CACNA2D1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CACNA2D1 CRISPRa kit - CRISPR gene activation of human calcium voltage-gated channel auxiliary subunit alpha2delta 1
Symbol:	CACNA2D1
UniProt ID:	P54289
Synonyms:	CACNA2; CACNL2A; CCHL2A; LINC01112; lncRNA-N3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)