

## Product datasheet for **GA100549**

### Human CACNB3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CACNB3 CRISPRa kit - CRISPR gene activation of human calcium voltage-gated channel auxiliary subunit beta 3
Symbol:	CACNB3
UniProt ID:	<a href="#">P54284</a>
Synonyms:	CAB3; CACNLB3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)