

Product datasheet for **GA100556**

Human S100G activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	S100G CRISPRa kit - CRISPR gene activation of human S100 calcium binding protein G
Symbol:	S100G
UniProt ID:	P29377
Synonyms:	CABP; CABP1; CABP9K; CALB3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)