

Product datasheet for **GA100563**

Human CALM3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CALM3 CRISPRa kit - CRISPR gene activation of human calmodulin 3
Symbol:	CALM3
UniProt ID:	P62158
Synonyms:	HEL-S-72; PHKD; PHKD3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)