

Product datasheet for **GA100570**

Human CAMK2D activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | CAMK2D CRISPRa kit - CRISPR gene activation of human calcium/calmodulin dependent protein kinase II delta |
| Symbol: | CAMK2D |
| UniProt ID: | Q13557 |
| Synonyms: | CAMKD |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)