

Product datasheet for **GA100603**

Human Caveolin 3 (CAV3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CAV3 CRISPRa kit - CRISPR gene activation of human caveolin 3
Symbol:	CAV3
UniProt ID:	P56539
Synonyms:	LGMD1C; LQT9; VIP-21; VIP21
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)