

Product datasheet for **GA100607**

Human CBFA2T3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CBFA2T3 CRISPRa kit - CRISPR gene activation of human CBFA2/RUNX1 partner transcriptional co-repressor 3
Symbol:	CBFA2T3
Components:	GA100607G1, CBFA2T3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA100607G2, CBFA2T3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA100607G3, CBFA2T3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q75081</u>
Synonyms:	ETO2; MTG16; MTGR2; ZMYND4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)