

Product datasheet for **GA100628**

Human Cyclin D3 (CCND3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CCND3 CRISPRa kit - CRISPR gene activation of human cyclin D3
Symbol:	CCND3
UniProt ID:	P30281
Synonyms:	cyclin D3; D3-type cyclin; G1/S-specific cyclin D3; OTTHUMP00000016390
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)