

Product datasheet for **GA100640**

Human CD1 (CD1D) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | CD1D CRISPRa kit - CRISPR gene activation of human CD1d molecule |
| Symbol: | CD1D |
| UniProt ID: | P15813 |
| Synonyms: | CD1A; R3; R3G1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)