

Product datasheet for **GA100659**

Human MS4A3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MS4A3 CRISPRa kit - CRISPR gene activation of human membrane spanning 4-domains A3
Symbol:	MS4A3
UniProt ID:	Q96HJ5
Synonyms:	CD20L; HTM4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)