

## Product datasheet for **GA100665**

### Human CD30 (TNFRSF8) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF8 CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 8
Symbol:	TNFRSF8
UniProt ID:	<u>P28908</u>
Synonyms:	CD30; D1S166E; Ki-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)