

Product datasheet for **GA100666**

Human CD153 (TNFSF8) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFSF8 CRISPRa kit - CRISPR gene activation of human TNF superfamily member 8
Symbol:	TNFSF8
UniProt ID:	P32971
Synonyms:	CD30L; CD30LG; CD153
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)